

AIR SCALERS FROM NSK  
**Ti-Max AIR SCALERS**

**NSK**  
**S900L/S900**

compact & powerful

NSK Ti-Max Air Scaler offers excellent power, low noise and connects to NSK standard air turbine coupling.



S900L Air Scaler handpiece delivers optic illumination direct to the clinical treatment site.



NSK Couplings

MODEL	ORDER CODE
<b>S900L</b>	T696

• Weight : 62g

**Optic**

<b>S900</b>	T695
-------------	------

• Weight : 65g

**Non-Optic**

• Includes 3 tips (S1, S2 & S3)  
 Tip Wrench and Tip Cover



KaVo®  
 MULTiflex LUX®

<b>S900KL</b>	T762
---------------	------

• Weight : 72g

**Optic**



Sirona®  
 Quick Coupling

<b>S900SL</b>	T774
---------------	------

• Weight : 68g

**Optic**




• Includes 3 tips (S1, S2 & S3)  
 Tip Wrench and Tip Cover







Powerful Partners®

# Air Scalars Tips for S900L / S900

## SCALING

MODEL	PRODUCTS	ORDER CODE
S1 Universal Type		Z252-411
S2 Sickle Type		Z252-413
S3 Perio Type		Z252-412

## PERIODONTICS

MODEL	PRODUCTS	ORDER CODE
S10 Universal Type Slim Line		Z252-420
S20 Root Planing Straight Tip		Z252-441
S21R Root Planing Right Curved Tip		Z252-442
S21L Root Planing Left Curved Tip		Z252-443

## RESTORATIVE (V-Tip)

Minimal Intervention	S75 V-Tip Holder		Z252-035
	V-G70 Diamond Tip		Y900-118
	V-G71 Diamond Tip		Y900-119
	V-G72 Diamond Tip		Y900-120
	V-G75 Diamond Tip		Y900-121
	V-G76 Diamond Tip		Y900-122
	Assorted Tip Kit	included 5 tips (V-G70, G71, G72, G75, G76)	Y900-158

Finishing / Trimming / Polishing	V-G77 Diamond Tip / Normal Grade		Y900-165
	V-G78 Diamond Tip / Fine Grade		Y900-166
	V-G79 Diamond Tip / Super Fine Grade		Y900-167
	Assorted Tip Kit	included 3 tips (V-G77, G78, G79)	Y900-168
Caries of Dentin	V-S1 Stainless Tip		Y900-161
	V-S3 Stainless Tip		Y900-162
	V-S33 Stainless Tip		Y900-163
	Assorted Tip Kit	included 3 tips (V-S1, V-S3, V-S33)	Y900-164


## Spare Parts



CR-2 Tip Replacement Wrench with Torque Limiter		Z221-075
Tip Holder		Z252-A080

Tip Cover L		Z217-850
Tip Cover S		Z217-851

## FlexiQuick Couplings for S900L / S900

**Titanium\***



	PTL-CL-4HV-T	P401-081
• LUX Coupling	• Optic Coupling for S900L	
	PTL-CL-M4-T	P434-081
• Midwest 4 hole	• Non-Optic Coupling for S900	

\*Only external body component is titanium

\*KaVo® and MULTiflex® are registered trademarks of Kaltenbach & Voigt GmbH & Co., Germany.  
Sirona® is a registered trademark of Sirona Dental Systems GmbH, Germany.



All NSK handpieces are built to have excellent durability against repeated autoclaving at temperatures of 135°C.



**NAKANISHI INC.**  
700 Shimohinata, Kanuma,  
Tochigi 322-8666, Japan  
TEL : +81 (0)289-64-3380  
FAX : +81 (0)289-62-5636  
www.nsk-inc.com

**TOKYO OFFICE**  
6F MY Building, Ueno 1-13-3,  
Taito-ku, Tokyo 110-0005, Japan  
TEL : +81 (0)3-3835-2892  
FAX : +81 (0)3-3835-2856



Westerbachstraße 58 D-60489 Frankfurt, Germany  
TEL : +49 (0)69 74 22 99 0 FAX : +49 (0)69 74 22 99 29  
www.nsk-europe.de



Dubai Airport Free Zone, Room 2W-310, 3F West Wing, Dubai, UAE  
TEL : +971-4-299-6644 FAX : +971-4-299-6664



Expect Perfection®