

MICRO SAW HANDPIECES FROM NSK
MICRO SAW HANDPIECES

NSK

MICRO SAW

sharp & effective



NSK Micro Saw Handpieces are designed to cut bones. Three types of actions are available - 1.8mm reciprocating action, 17° oscillating action and 3° sagittal ation, to perform easy and effective bone cutting operations in oral surgery. The NSK Micro Saw Handpiece range is designed to be ultra thin to minimize damage of bone during cutting operations, and to allow excellent treatment site access and visibility.



Powerful Partners®




Reciprocating 1.8 mm vertical reciprocating action



MODEL	ORDER CODE
SGR-E	H173

- Maximum Speed Allowance : 30,000 min⁻¹
- With External Spray Nozzle



Oscillating 17° horizontal reciprocating action



MODEL	ORDER CODE
SGO-E	H174

- Maximum Speed Allowance : 15,000 min⁻¹
- With External Spray Nozzle



Sagittal 3° horizontal reciprocating action






MODEL	ORDER CODE
SGT-E	H175







- Maximum Speed Allowance : 20,000 min⁻¹
- With External Spray Nozzle

*Do not exceed the Maximum Speed Allowance.

SGR-E Blades for Reciprocating

PRODUCTS	MODEL	ORDER CODE
 Blade thickness 0.35 mm	SGR-1 Pack of 10 blades	Y900-072
 Blade thickness 0.35 mm	SGR-2 Pack of 10 blades	Y900-073
 Blade thickness 0.35 mm	SGR-3 Pack of 10 blades	Y900-074

SGO-E Blades for Oscillating

PRODUCTS	MODEL	ORDER CODE
 Blade thickness 0.3 mm	SGO-1 Single blade	H174-034
 Blade thickness 0.3 mm	SGO-2 Single blade	H174-044
 Blade thickness 0.3 mm	SGO-3 Single blade	H174-031
 Blade thickness 0.3 mm	SGO-4 Single blade	H174-041
 Blade thickness 0.35 mm	SGO-5 Single blade	H174-032
 Blade thickness 0.35 mm	SGO-6 Single blade	H174-042

SGT-E Blades for Sagittal

 Blade thickness 0.35 mm	SGT-1 Pack of 10 blades	Y900-075
 Blade thickness 0.35 mm	SGT-2 Single blade	Y900-076

Specifications are subject to change without notice.



All NSK handpieces are autoclavable up to 135°C.



www.nsk-inc.com

NAKANISHI INC.
700 Shimohinata
Kanuma, Tochigi
322-8666, Japan
TEL:+81 (0)289-64-3380
FAX:+81 (0)289-62-5636

NSK TOKYO Office
6F MY Building
Ueno 1-13-3, Taito-ku
Tokyo 110-0005, Japan
TEL:+81 (0)3-3335-2892
FAX:+81 (0)3-3335-2856

NSK EUROPE GmbH
Westerbachstraße 58
D-60489, Frankfurt
Germany
TEL:+49 (0)69 74 22 99 0
FAX:+49 (0)69 74 22 99 29

NSK FRANCE SAS
71 rue de Suresnes,
92380 Garches
France
TEL:+33 (0)1-4795-4517
FAX:+33 (0)1-4741-3516

NSK MIDDLE EAST
Dubai Airport Free Zone,
Room 2W-310,
3F West Wing, Dubai, UAE
TEL:+971-4-298-6644
FAX:+971-4-299-6664



Expect Perfection®