

NSK

OPTIC AIR MOTOR

Ti205L

smooth and quiet



Powerful Partners®

Ti205L Optic Air Motor

ISO E-Type

Optic



135°C
SSS

The optic Ti205L air motor delivers powerful, smooth and quiet rotation with loads of torque at all speeds. Speed can be varied to a maximum of 20,500min⁻¹ in the forward direction. Ti205L air motor is highly durable and fully autoclavable up to 135°C. Ti205L connects directly to all brands of E-type optic and non-optic low speed handpieces.



360° Rotation with Illumination

Ti205L allows full 360° swivel of low speed handpieces, without the loss of illumination. The NSK anti-twist system prevents hose drag and minimises hand fatigue.



MODEL

ORDER CODE

Ti205L

M190

- 1:1 Direct Drive
- Max 20,500min⁻¹
- Internal Spray
- Lux Coupling
- NSK non-retraction valve

Ti-Max TITANIUM Ti Series Popular Attachments

Optic



135°C
SSS

NSK now uses solid titanium extensively for the latest low speed handpiece models. Ti-Max series low speed handpieces incorporate a new "slim-line" neck design that allows clear operator vision direct to the operating field while the integrated cellular optics provides clear illumination. Connecting NSK Ti-Max low speed handpieces to the Ti205L optic air motor creates the ideal combination for all low speed procedures. Other brands of optic and non-optic low speed handpieces can also be connected to the NSK Ti205L air motor.



Ti 95L 1:5 Speed Increasing

MODEL

ORDER CODE

Ti 95L

C448

- 1:5 Increasing
- Quattro spray
- For FG burs (ø2.35)



Ti 15L 4:1 Speed Reduction

MODEL

ORDER CODE

Ti 15L

C444

- 4:1 reduction
- Single spray
- For CA burs (ø2.35)



Ti 25L 1:1 Drive

Ti 25L

C446

- 1:1 Drive
- Single spray
- For CA burs (ø2.35)



Ti 65L Straight

Ti 65L

H248

- 1:1 Drive
- Single spray
- For HP burs (ø2.35)

NAKANISHI INC.
www.nsk-inc.com
700 Shimohinata Kanuma-shi
Tochigi 322-8666,
Japan

TEL : +81 (0)289-64-3380
FAX : +81 (0)289-62-5636

NSK Europe GmbH
www.nsk-europe.de
Westerbachstraße 58
D-60489 Frankfurt,
Germany

TEL : +49 (0)69 74 22 99 0
FAX : +49 (0)69 74 22 99 29

NSK United Kingdom Ltd.
www.nsk-inc.com
Gateway1000, Arlington Business Park,
Whittle Way, Stevenage,
Hertfordshire, SG1 2FP, UK

TEL : +44 800 634 1909
FAX : +44 800 634 1910

NSK Oceania Pty Ltd.
www.nsk-inc.com
Unit 22, 198-222 Young St.
Waterloo, Sydney,
NSW 2017 Australia

TEL : +61 (0)2 8306 3000
FAX : +61 (0)2 8306 3030

NSK America - Canadian Branch
www.nsk-inc.com
1054 Centre St. - Suite 283
Thornhill, Ontario, Canada
L4J 8E5

TEL/FAX : +1 905 780 6757
Toll Free : +1 866 780 6757

NSK Middle East
www.nsk-inc.com
Dubai Airport Free Zone,
Room 2W-310,
3F West Wing, Dubai, UAE

TEL : +971 (0)4-299-6644
FAX : +971 (0)4-299-6644