

Ti-Max M40 FROM NSK
ELECTRIC MICROMOTOR

NSK
M40

powerful, reliable and flexible

M40 Easily Fits on Bien-Air®
4VLM, 4VR 400 and B4VLR



Powerful Partners®

Ti-Max M40

Optic

ISO E-Type Internal Water Spray

Connect Direct to Bien-Air® 4VLM, 4VR 400 and B4VLR hoses.

The new NSK M40 micromotor employs the advanced technology of coreless DC motor. This technology smoothly delivers powerful torque and also minimizes any heat generation. M40 provides a small, yet high performance and quiet operation. M40 is able to connect to Bien-Air® short rotary connector. The clinician can benefit from a wide selection of instruments.




Features


- The ideal combination of speed and torque
- Minimum heat generation
- Quiet smooth operation from 60 to 40,000min⁻¹
- Optic micromotor
- Lightweight



NSK Ti-Max M40




NSK



Bien-Air®

4VLM / 4VR 400 / B4VLR
(Made by Bien-Air Dental S.A.)



Ti-Max M40

E317-051

- 60~40,000min⁻¹
- M40 micromotor only
- Optic
- Includes replacement TM bulbs and O-rings

Specifications

- Speed : 60~40,000min⁻¹
- Dimensions : L103 x ø21 (mm)
- Weight : 106g

Bien-Air® are registered trademarks of Bien-Air Dental S.A., Switzerland.



www.nsk-inc.com

NAKANISHI INC
700 Shimohinata
Kanuma, Tochigi
322-8666, Japan
TEL:+81 (0)289-64-3380
FAX:+81 (0)289-62-5636

NSK TOKYO OFFICE
6F MY Building
Ueno 1-13-3, Taito-ku
Tokyo 110-0005, Japan
TEL:+81 (0)3-3835-2892
FAX:+81 (0)3-3835-2856

NSK EUROPE GmbH EC REP
Westerbachstraße 58
D-60489, Frankfurt,
Germany
TEL:+49 (0)69 74 22 99 0
FAX:+49 (0)69 74 22 99 29

NSK FRANCE SAS
19 avenue de Villiers
75017 Paris,
France
TEL:+33 01 56 79 59 80
FAX:+33 01 56 79 59 81

NSK MIDDLE EAST
Dubai Airport Free Zone,
Room 2W-310,
3F West Wing, Dubai, UAE
TEL:+971-4-299-6644
FAX:+971-4-299-6664

Specifications are subject to change without notice.

REF: PR-D329E Ver.1 06.04.05.Y